



LED Headlight Systems

BR Surgical medical headlights are engineered for long-lasting service, unbeatable performance and superior comfort.



BR900-3191

Wireless LED Headlight with 2 rechargeable batteries, single bay charger and carrying bag

- Integrated LED light source and deluxe headband with crown and cranial support
- Lightweight and cable-free
- Color Temperature: 4500K
- Color Rendition Index (CRI): 75
- Typical LED Life: >50,000 hours
- Variable Module Spot: 20mm-110mm at 16" (40cm)
- Side Mounted Variable Intensity Control: 0-30,000 Lux measured at 16" (40cm)
- Multidirectional light adjustment
- Battery Life: 4 hours continuous operating time
- Battery Recharge Time: 4.5 hours
- Warranty: Three (3) years on LED Headlight and one (1) year on Batteries

SYSTEM INCLUDES

- 1 LED headlight with rear cranial support headband
- Side mounted battery holster and attenuator
- 2 Rechargeable batteries
- 1 Single bay charger
- 1 Padded carrying bag with shoulder strap



BR900-3192

LED Headlight with battery pack, 2 rechargeable batteries, single bay charger and carrying bag

- Integrated LED light source and deluxe headband
- Lightweight
- Color Temperature: 4500K
- Color Rendition Index (CRI): 75
- Typical LED Life: >50,000 hours
- Variable Module Spot: 20mm-110mm at 16" (40cm)
- Variable Intensity Control: 0-70,000 Lux measured at 16" (40cm)
- Headband with LED Module Weight: 9 oz. (225g)
- Battery and Holster Weight: 12 oz. (340g)
- Battery Life: 2 hours 45 minutes at full intensity
5 hours 30 minutes at half intensity
- Battery Recharge Time: 3.5 hours
- Battery charger accepts 100-240V AC, 50/60 Hz input for worldwide operation
- Warranty: Three (3) years on LED Headlight and one (1) year on Batteries

SYSTEM INCLUDES

- 1 LED headlight with deluxe headband
- 1 battery holster
- 2 Rechargeable batteries
- 1 Single bay charger
- 1 Padded carrying bag with shoulder strap