

TOE NAIL PROCEDURE KIT

BR980-74000



BR Surgical provides a **Sterile-Disposable Toe Nail Procedure Kit** for in-grown toe nail removal and other nail removal procedures. The tray includes (5) disposable instruments necessary to perform the procedure.

Simple, sterile, and ready to use on your everyday toe nail procedures.

- 1) Kelly Forceps (straight):** Removing/pulling back tissue around the nail
- 2) Adson Forceps 1x2:** Grasping tissue or the nail
- 3) Iris Scissors (straight):** Cutting surrounding nail tissue
- 4) Probe with Eye:** Elevating/manipulating and separating the nail
- 5) Nail Nipper:** Cutting and trimming the nail

NOTE: This **BR980-74000** sterile kit is sold with (5) trays per box.

CPT Codes:

- 11730** - Avulsion of nail plate, partial or complete, simple; single
- 11732** - Avulsion of nail plate, partial or complete, simple; each additional nail plate (list separately in addition to code for primary procedure)
- 11750** - Excision of nail and nail matrix, partial or complete (e.g., ingrown or deformed nail), for permanent removal
- 11765** - Wedge excision of skin of nail fold (e.g., for ingrown toenail)



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com
Office: 904.642.1366 • Toll Free: 888.642.1366 • Fax: 904.642.1368 • www.BRSurgical.com