



ENT Treatment Cabinet

BR900-7951

BR Surgical offers a compact and affordable ENT Treatment Cabinet

- Large, easy to clean laminate work surface
- Drawers with organized inserts
- Workstation ON/OFF switch for suction and pressure pump control
- Pressure and vacuum hose holder on side of cabinet
- 1200cc suction canister
- Premium suction and pressure pumps

INCLUDES:

- (1) Atomizer Bottle
- (1) Nebulizer Bottle
- (2) Stainless Steel Jars with Lid

DIMENSIONS:

- Overall Height: 47 1/2"
- Height to Work Surface: 33 1/4"
- Width: 24 1/2"
- Depth: 20 1/4"



Ask us about our line of ENT Chairs, Treatment Cabinets, Microscopes, Headlights, and Stroboscopes.

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com